Ciclosporin A, Tacrolimus, Sirolimus und Everolimus:

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Literatur:

- (1) Am J Kidney Dis. 1996 Aug;28(2):159-72. Immunosuppressants: cellular and molecular mechanisms of action. Suthanthiran M, Morris RE, Strom TB.
- (1b) Clinical Chemistry. 2002;48:2225-2231. Delayed Cytokine mRNA Expression Kinetics after T-Lymphocyte Costimulation: A Quantitative Measure of the Efficacy of Cyclosporin A-based Immunosuppression. Christoph Härtel1, Lutz Fricke2, Nina Schumacher1, Holger Kirchner1 and Michael Müller-Steinhardt1
- (2) Science. 1989 Dec 22;246(4937):1617-20. Cyclosporin A specifically inhibits function of nuclear proteins involved in T cell activation. Emmel EA, Verweij CL, Durand DB, Higgins KM, Lacy E, Crabtree GR.
- (2b) J Lab Med 2003;27 (5/6):222-227. Drug Monitoring of Sirolimus and Everolimus. Victor W. Armstrong, Frank Streit
- (3) Cancer Biol Ther. 2003 May-Jun;2(3):222-32. Rapamycin: mechanism of action and cellular resistance. Huang S, Bjornsti MA, Houghton PJ
- (4) Clin Breast Cancer. 2003 Jun;4(2):126-37. Mammalian target of rapamycin: a new molecular target for breast cancer. Mita MM, Mita A, Rowinsky EK
- (5) Mol Immunol. 2003 Jul;39(17-18):1073-7. The mosaic of immunosuppressive drugs. Masri MA
- (8) EMBO J. 1995 Jun 15;14(12):2771-83. Targets of immunophilin-immunosuppressant complexes are distinct highly conserved regions of calcineurin A. Cardenas ME, Muir RS, Breuder T, Heitman J

Control of Assays for Immunosuppressive Drugs

Dr. David W Holt, Director, Analytical Unit, St George's Hospital School, London SW17 0RE
Contact: d.holt@sgghms.ac.uk


The immunosuppressive drugs are characterised by substantial between-patient differences in their pharmacokinetics and a broad spectrum of unwanted effects. Cyclosporin and tacrolimus, still the mainstay of drug therapy for the avoidance of graft rejection, are routinely measured in blood to assist in their optimal prescription. Recently, attention has focused on a newer drug, sirolimus, for which drug monitoring in blood is a regulatory requirement in Europe.

There are no established reference methods for the measurement of the immunosuppressive drugs, nor are there certified reference materials for calibration. As a result, between-method assay performance varies substantially. This is due both to errors in calibration, and to differences in specificity of the antibodies used in the immunoassays. Data on method performance can be judged by external proficiency testing material, as presented on our web site (www.bioanalytics.co.uk), and by reference to consensus documents on immunosuppressive drug monitoring. A joint Working Group sponsored by the International Federation of Clinical Chemistry (IFCC) and the International Association of Therapeutic Drug Monitoring and Clinical Toxicology is reviewing current methods and monitoring strategies (1).

There are some interesting developments in immunosuppressive drug therapy, related to both the introduction of new drugs and the use of new drug combinations. These developments pose analytical challenges for the laboratory related to sensitivity

and selectivity. A driving force for innovation in methodology has been the introduction of sirolimus since, currently, HPLC is the only technique available for its measurement. More than 100 centres world-wide are now measuring the drug by HPLC, about half of these using mass-spectrometric detection (LC/MS). This technique

has also underpinned all the clinical studies of the sirolimus analogue, everolimus, to be introduced in Europe early in 2004, and is being used by an increasing number of centres for the measurement of cyclosporin and tacrolimus.

and selectivity. A driving force for innovation in methodology has been the introduction of sirolimus since, currently, HPLC is the only technique available for its measurement. More than 100 centres world-wide are now measuring the drug by HPLC, about half of these using mass-spectrometric detection (LC/MS). This technique

The need for centres to prepare in-house calibrators for use in LC/MS assays has focused attention on the need for well-defined control material. This is particularly so for sirolimus, since there are no immunoassay results with which to compare locally prepared calibration material. Recently, we have tested a set of four freeze-

dried whole blood controls with defined-values for the measurement of sirolimus, using LC/MS (Chromsystems GmbH, München, Germany). The results are illustrated below and show good agreement between the mean of 5 replicates of each control, analysed in this laboratory, and the mean defined value. The average deviation from the expected value was 6.5%. A blank matrix and calibrator are also available from the same source. All these materials also contain defined values of cyclosporin and tacrolimus.

Even for immunoassays there is interest by manufacturers in defining their own calibrators against well sourced control material. Currently, we are supplying services to diagnostics manufacturers for the assessment of calibrator material using validated LC/MS assays, and by supplying proficiency testing material previously assayed on multiple sites.

Thus, in a continuously evolving field, there is a need for constant vigilance to ensure that assays for these critical dose drugs are both accurate and precise.

Literature:

- (1) Ther Drug Monit, 24, 59-67, 2002. International Federation of Clinical Chemistry/International Association of Therapeutic Drug Monitoring and Clinical Toxicology Working Group on immunosuppressive drug monitoring. Holt DW, Armstrong VW, Griesmacher A, Morris RG, Napoli KL, Shaw LM.

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Literatur

- (1) Drug Safety, 1990, 5, 359-383. An introduction to the clinical toxicology of volatile substances, R. J. Flanagan, M. Ruprah, T. J. Meredith, J. D. Ramsey,
- (2) Medical Hypotheses, 2000, 54, 619-623. New aspects in genotoxic risk assessment of styrene exposure - a working hypothesis, B. Marczynski, M. Peel, X. Baur
- (3) Pathology, 1994, 26, 301-309. Xylene: its toxicity, measurement of exposure levels, absorption, metabolism and clearance, J. L. Langman
- (4) Toxicology Letters 1995, 77, 85-91. Exposure to complex mixtures: implications for biological monitoring, M. Ikeda

Determination of solvent Metabolites in Urine

Dr. Claudio Minoia,
Laboratory for Environmental
and Toxicological Testing,
"Salvatore Maugeri" Foundation, Pavia, Italy

Aromatic solvents, such as toluene, xylene and styrene, are widely used in several applications, including industrial solvents and chemical intermediates. Since these compounds are often used in combination with each other, workers are usually exposed to a mixture of solvents, which can also be due to the presence of significant amounts of contaminants. Occupational exposure is commonly evaluated by biological monitoring with simultaneous measurement of the working environment. Since solvents are mainly absorbed by inhalation, air concentrations of these compounds are measured, and urinary metabolites are also determined. The results from both environmental and biological monitoring are then compared for example with the corresponding Threshold Limit Values (TLVs) and Biological Exposure Indexes (BEIs) developed by the American Conference of Governmental Industrial Hygienists (ACGIH).

Nevertheless, detectable concentrations of some of these metabolites, such as hippuric acid, can be found in the urine of non-exposed subjects. The Reference Values in the general population should be therefore taken into account, and the influence of different variables (diet, lifestyle, etc.) on the urinary levels of these substances should be evaluated.

The main end products of styrene metabolism are urinary mandelic acid (MA) and phenylglyoxylic acid (PGA). In the occupational studies, end-of-shift and next-morning MA + PGA are used as indicators of chronic styrene exposure. According to the ACGIH, the air level of styrene should not exceed 20 ppm (86 mg/m³). This value is strictly related to the next-morning urinary mandelic acid and

phenylglyoxylic acid, which should be less than 300 mg/g creatinine and 100 mg/g creatinine, respectively. As far as end-of-shift urinary levels are concerned, the mandelic and phenylglyoxylic acid concentrations should not be higher than 800 mg/g creatinine and 240 mg/g creatinine, respectively. In contrast, in the general population, average levels of 3.8 mg mandelic acid/g creatinine and 4.3 mg phenylglyoxylic/g creatinine are measured.

Toluene is commonly used as an industrial solvent for the manufacturing of paints, chemicals, pharmaceuticals, and rubber and a TLV-TWA of 50 ppm (188 mg/m³) has been established by the ACGIH. Toluene exposure levels can be determined from urinary hippuric acid levels, which is its major urinary metabolite. Nevertheless, when workers are exposed to low levels of airborne toluene (short exposures or TWA concentrations below 10 ppm), correlation coefficients between end-of-shift levels of hippuric acid and toluene concentrations in

air, are not statistically significant. Urinary hippuric acid levels are in this case very close to those measured in the general population (average levels 300-400 mg/g creatinine, ranging from 50 to 2000 mg/g creatinine). This means that there are a number of hippuric acid precursors in the environment, such as benzoic acid in food (e.g. plums, preserved food containing fish or eggs), which makes diet one of the major confounding factors.

Xylene is a colourless liquid that catches on fire easily. Commercial or mixed xylenes are composed of three isomers (meta-, ortho- and para-xylene) and are extensively used in the chemical industry. Occupational exposure to xylene can be assessed by determining urinary concentrations of its metabolites, m-, o- and p-methylhippuric acids. The ACGIH has established a TLV-TWA for airborne xylene of 100 ppm (350 mg/m³), whereas end-of-shift level of the three acids should be less than 1500 mg/g creatinine. The reference values in

the general population are 0.59, 18.9, and 0.43 mg/g creatinine for meta-, ortho-, and para-methylhippuric acid, respectively. Since, as mentioned above, both general population and workers are usually exposed to mixtures of solvents, a simple and sensitive method for the quantification of even low urinary levels of their major metabolites is needed.

A reagent kit for the simultaneous determination of mandelic, hippuric, phenylglyoxylic, ortho-, meta-, para-methylhippuric acids has been developed by Chromsystems – Diagnostic by HPLC. The sample preparation is very simple and just 10 ml of urine sample are required. A chromatographic run of 18 minutes allows the separation of six metabolites which can be quantified by the inclusion of the internal standard. The results can be then compared with the BEIs established by the ACGIH.

This reagent kit is a time- and cost-saving tool for risk assessment related to the exposure to toluene, xylene, and styrene.

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Literatur:

- (1) Curr Med Chem. 2004 Feb;11(3):343-58. New antipsychotics and schizophrenia: a review on efficacy and side effects. Serretti A, De Ronchi D, Lorenzi C, Berardi D.
- (2) Curr Med Chem. 2004 Feb;11(3):279-96. Atypical antipsychotics: pharmacokinetics, therapeutic drug monitoring and pharmacological interactions. Raggi MA, Mandrioli R, Sabbioni C, Pucci V.
- (3) J Clin Psychiatry. 2004 Sep;65 Suppl 9:3-8. Historical perspective on movement disorders. Friedman JH.
- (4) www.klein-kreienkamp-oehler.de/FG_Psychiatrie.htm, Neuroleptika, ihre Indikation, Wirkung, Risiken und Nebenwirkungen, Wirkungsweise
- (5) JAMA. 2003 Nov 26;290(20):2693-702. Effectiveness and cost of olanzapine and haloperidol in the treatment of schizophrenia: a randomized controlled trial. Rosenheck R, Perlick D, Bingham S, Liu-Mares W, Collins J, Warren S, Leslie D, Allan E, Campbell EC, Caroff S, Corwin J, Davis L, Douyon R, Dunn L, Evans D, Frecka E, Grabowski J, Graeber D, Herz L, Kwon K, Lawson W, Mena F, Sheikh J, Smelson D, Smith-Gamble V.
- (6) Pro Mente Sana. 2003 Jan, 16-18. Atypische Neuroleptika: Unwahrheiten für schlichte Gemüter. Lehmann P. www.antipsychiatrie.berlinet.de/artikel/gesundheit/atypische.htm
- (7) Int Clin Psychopharmacol. 2004 Jul;19(4):251-253. Weight gain during long-term treatment with olanzapine: a case series. Haberer EM, Rittmannsberger H.
- (8) Can J Psychiatry. 2004 May;49(5):343. Tardive dyskinesia associated with olanzapine in a neuroleptic-naïve patient with schizophrenia. Bhanji NH, Margolese HC.
- (9) Pharmacopsychiatry. 2004 May;37(3):131-4. Olanzapine- and clozapine-induced stuttering. A case series. Bar KJ, Hager F, Sauer H.
- (10) World J Biol Psychiatry. 2004 Apr;5(2):73-82. Atypical antipsychotics and diabetes mellitus. Schwenkreis P, Assion HJ.

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Instruments & Chemicals GmbH